



# Sample Preparation

## Module 1: Overview

# Sequencing Workflow



**Sample  
Preparation**



**Cluster  
Generation**



**Sequencing**



**Data Analysis**




# Library Preparation


DNA  
(0.1-5.0 µg)



Library  
preparation



Cluster growth



Sequencing


# Library Preparation



**Sample  
Preparation**

- ▶ Prepares sample nucleic acid for sequencing
  - Fragmenting
  - Generates double-stranded DNA (if necessary)
  - Flanks with Illumina adapters
- ▶ All preparation ends with the same general template structure
  - Double-stranded DNA flanked by adapters
  - Variables include
    - Insert Size
    - Adaptor type
    - Index

# Library Preparation



## Library preparation:


Attach adapters to both ends of fragmented nucleic acid



BROAD  
INSTITUTE

illumina®

# Genomic DNA Protocol Overview



- ▶ Most popular protocol
- ▶ Basis for all other protocols

# Starting Material Input


- ▶ Starting Materials

- Quantify
- Quality Control

- ▶ Methods for Quantification:



NanoDrop\*



PicoGreen® \*





Qubit™

- ▶ Methods for Quality Control

- Visualization on gel
- Agilent Bioanalyzer\* trace especially useful for RNA (RIN)

\*Methods utilized at the Broad

# Step 1: Fragment DNA





# Step 1: Fragment DNA

Purified Genomic DNA

Fragment DNA

↓ Fragments < 800 bp

Repair Ends

↓ Blunt End Fragments with 5' –  
Phosphorylated ends

Add an “A” to the 3' Ends

↓ 3' – dA Overhang

Ligate Paired-End  
Adapters

↓ Adapter-Modified Ends

Size Select on Gel



↓ 300 – 600 bp Fragments

PCR

↓ Amplified DNA with Adapters


QC Library

Genomic DNA Library



- ▶ Nebulizer
  - Very inexpensive
  - Works well for 1-5  $\mu$ g starting input
- ▶ Other methods available
  - Covaris™
  - Sonication
  - HydroShear®
  - Enzymatic
  - Others



# Step 1: Fragment DNA




	Nebulizer	Covaris™	HydroShear®	Sonication	Chemical	Enzymatic
DNA	✓	✓	✓	✓		
mRNA					✓	
DGE						✓
small RNA						
CHIP-Seq		✓		✓		
Mate Pair	✓		✓			

- Fragmentation methods have benefits and limitations
- Can be application-specific
- Depends on desired end-product

# Step 2: Repair Ends





# Step 2: Repair Ends




- ▶ Converts overhangs
  - Blunt ends
  - Phosphorylates 5' End
- ▶ Reagents:
  - dNTP, T4 DNA pol, Klenow
  - Kinase/ATP (T4 PNK)
- ▶ Simple enzymatic reaction
- ▶ 30 minutes incubation

# Step 3: Add “A” Overhang





# Step 3: Add “A” Overhang




- ▶ Adds
  - ‘A’ base to the 3' end of the blunt phosphorylated DNA fragments
- ▶ Prevents
  - Formation of adapters dimers
  - Concatemers
- ▶ Reagents
  - 1 mM dATP, Klenow exo (3' to 5' exo minus)

# Step 4: Ligate Adapters





# Step 4: Ligate Adapters



- ▶ Same library structure regardless of Application
  - gDNA
  - RNA
  - miRNA
  - ChIP-seq
  - Exon Pull-down
  - etc.



# Step 4: Ligate Adapters




- ▶ 10:1 molar ratio of adapter to genomic DNA insert
  - Assumes 5 µg initial input
  - Titrate for different insert amounts

# Choosing a Library Type

## ▶ Single read libraries:

- Unidirectional Sequencing
- Single Read Flowcells ONLY
- Counting applications: ChIP or low coverage resequencing projects




## ▶ Paired end libraries:

- Uni- OR Bi-directional (paired reads)
- Paired End Flowcells; Single: Unidirectional only
- Most applications, #1 whole genome shotgun assembly
- Tailor insert size and distribution per project:
  - Tight size distribution – Assembly, structural rearrangement detection
  - Wide distribution libraries - Resequencing, high coverage



## ▶ Multiplex Paired End (aka Indexing or Barcoding)

- Uni- OR Bi-directional
- Allows multiple libraries per lane
- 12 Index tags available x 8 lanes = 96 libraries per flowcell




# Flowcell Compatibility

## Single Read Flowcell

Periodate Linearization:

- Cleaves off most of the P5 oligo



5' P-GATCGGAAGAGCTCGTATGCCGTCTTCTGCTTG  
5' ACACTCTTTCCCTACACGACGCTCTTCCGATCT


SR Adaptors

## Paired End Flowcell

Linearization 1 Enzyme:

- Leaves part of P5 oligo intact for Read 2 resynthesis



Linearization 2 Enzyme




5' P-GATCGGAAGAGCGGTTCAGCAGGAATGCCGAG  
5' ACACTCTTTCCCTACACGACGCTCTTCCGATCT

PE Adaptors

# Step 5: Gel Purify Ligation Products




- ▶ 1µg of DNA Starting Material
  - 300/600bp area excised




- ▶ 5µg of DNA Starting Material
  - 300/600bp area excised

# Step 6: PCR Enrich Ligation Products




- ▶ Selectively enrich DNA fragments with adapter molecules on both ends
- ▶ Adds additional sequences to the end of the adapters for hybridization
- ▶ Amplifies the amount of DNA in the library
- ▶ Reagents:
  - Adapter Specific Primers
  - High Fidelity Enzyme (Phusion)


# Step 7: QC Library



- ▶ Quantitate by NanoDrop or Qubit
  - NanoDrop can overestimate concentration
  - Qubit or PicoGreen gives more accurate measurement
- ▶ QC by Agilent Bioanalyzer®
  - Gives size confirmation
  - Visualizes unwanted products





# Step 7: QC Library



## ► Agilent Bioanalyzer Results

- Lane 1: 300 bp library
- Lane 2: 600 bp library





# Other DNA Applications

- ▶ Multiplexing (aka Indexing/Barcoding)
- ▶ Large Insert Mate Pair




# Multiplex Paired End Libraries

Index Tag is incorporated during Sample Prep PCR, using a 3-primer system



# Multiplex Paired End Libraries

Multiplex Sequencing Utilizes 3 Sequencing Reads



# Long Insert Mate Pair

- ▶ Allows gap sizes of 2-5kb
- ▶ Useful for whole genome shotgun assembly; connecting contigs
- ▶ Detection of structural/chromosomal rearrangements





Figure 1 Mate Pair Library Preparation Overview

# RNA Applications Sample Prep

- ▶ mRNA-Seq
- ▶ Small RNA Seq


# mRNA-Seq Sample Prep

- ▶ Isolate poly-A containing mRNA
  - Magnetic oligo dT beads
- ▶ Randomly fragment RNA
- ▶ Random prime mRNA → cDNA
- ▶ Make 2nd strand cDNA
- ▶ Ligate sequencing adapters
- ▶ Size select on gel
- ▶ Enrich with 15 cycles of PCR




# Diagram of Ligation Reaction

- ▶ Total RNA
- ▶ 3' adapter ligation
- ▶ 5' adapter ligation
- ▶ RT-PCR
- ▶ Sequencing



# Small RNA V1.5 Sample Prep

- ▶ 6 hour library preparation
  - Shortest time from library to sequence
  - 2.5 hrs hands on time (6 hr total)
  - Single gel excision
- ▶ NO need for small RNA purification!
  - Starts directly from total RNA
- ▶ Reduced RNA input requirements
  - 1.0 µg total RNA input
    - Titrated to 200 ng total
- ▶ Strand specific information
  - Discover and profile microRNA
  - Sequence other non-coding RNAs
  - Can be used with stranded mRNA-Seq



# Sample Preparation Best Practices



## Sample Preparation

- ▶ Be familiar with the protocol before starting it
  - Have all user-supplied reagents and equipment
- ▶ Use high-quality starting material
  - QC and quantify material before starting protocol
  - Aim for middle to high end of starting material requirements if possible
  - Very important for RNA applications
- ▶ Use proper material handling techniques
  - Prevent contamination and degradation
  - Particularly important for RNA applications



# Sample Preparation Best Practices



**Sample  
Preparation**

- ▶ Use positive controls if necessary for troubleshooting
- ▶ Quantitate and QC final library
  - This is essential for cluster lawn quality and sequencing success



GA BOOT CAMP


# Sample Preparation

Best Practices from The Broad Institute

# Module Objectives – Sample Preparation


- ▶ **Overview of Broad's sample preparation workflow**
- ▶ **How we track samples through the lab**
- ▶ **Broad developed modifications to sample preparation and their application in our protocols:**
  - ▶ Shearing methods
  - ▶ Enzymatic reaction clean up methods
  - ▶ Size selection methods
    - ▶ Both agarose gel and gel-free methods
  - ▶ PCR enrichment protocols
- ▶ **Quality control methodology for sample preparation**
  - ▶ Overview of sequencing metrics that can be used to assess sample preparation quality

# Sample Prep Workflow at the Broad



# Sample Prep Workflow at the Broad


*All genomic DNA tracked and aliquoted by our Sample Repository group*



# Sample Tracking

# Sample Tracking

*All genomic DNA tracked and aliquoted by our Sample Repository group*



# Sample Tracking in a High Throughput World

- ▶ DNA stored in 2D barcoded tubes
  - ▶ Source Genomic DNA
  - ▶ Unenriched libraries
  - ▶ Enriched libraries
  - ▶ Normalized libraries
- ▶ All transfers and details are tracked via web-based LIMS
- ▶ Tube picking system error-proofs selecting libraries



**BROAD INSTITUTE** Squid Index User: Andrew Barry  
[Logout](#) [Email](#) [Help](#)

**Projects** **Clones** **L.C** **Pipeline** **Detection** **Operations** **Admin** **Finishing** **Lims** **454** **Sample Tracking** **Solexa** **HT PCR**

---

**1 Source Receptacle(s):** 0061958247 Transferred Volume (uL): 4

**Barcode:** 0061958247  
**Material Type:** Solexa Library -> Normalized Library

**Content Details**  
Name: Solexa-1799  
Type: Solexa Library  
Remaining Volume: 27.7 uL  
Concentration: 0.9995 ng/uL  
Receptacle Barcode: 0061958247  
Current Buffer: In Sequencing  
Created On: 09-Oct-2007  
Created By: alexall  
Status: Active  
Description: Aliquot for Production Sequence  
GSSR ID: 17309.0  
Organism: Human  
SOP: Broad Solexa LC Protocol v2.0.doc  
Key Sequence: —  
Adaptor A: —  
Adaptor B: —  
Reference Sequence: Homo\_sapiens\_assembly18  
Library Lot #: —  
Fragment Size: 150.0  
Library Type: Normalized Library  
Analysis Type: Standard  
Absorbion 260: 0.0  
Absorbion 280: 0.0  
Concentration (nM): 10.149  
Single Stranded: No


**Transfer History:**  
**Backwards History**  
17309.0  
↓  
0061958234  
↓  
0061958247  
**Forwards History**  
0061958247  
↓  
9736\_Lane\_1

**1 Destination Receptacle(s):** 000000000 [Hide Details](#) [remove](#)

**Material Type:** Solexa Library -> Normalized Library  
**Volume (uL):** 4  
**Concentration (ng/uL):** 1.2  
**Reference Sequence:** Homo\_sapiens\_assembly18  
**Fragment Size (Insert + Adaptor):** 150  
**Library Kit Lot #:** 45678  
**Analysis Type:** Standard  
**SOP:** Broad Solexa LC Protocol v2.0.doc  
**Description:** Aliquot for Production Sequence



# Sample Tracking at Lower Throughputs




	A	B	C	D	E	F	G	H	I	J
1	<b>Library Tracking Sheet - October 2007</b>									
2										
3	<b>Project ID</b>	<b>Technician</b>	<b>Date LC Completed</b>	<b>Unenriched Library</b>	<b>Unenriched Location</b>	<b>Enriched Library</b>	<b>Enriched Location</b>	<b>Normalized Library</b>	<b>Normalized Location</b>	<b>Loading pM from QPCR</b>
4	E. coli MG1655 WGS	Holly	10/8/2007	Solexa-15901	Unrchd Box #3	Solexa-15912	Enrchd Box #7	Solexa-15923	Norm Box #11	3.63
5	E. coli MG1655 WGS	Holly	10/8/2007	Solexa-15902	Unrchd Box #3	Solexa-15913	Enrchd Box #7	Solexa-15924	Norm Box #11	0.50
6	E. coli MG1655 WGS	Holly	10/8/2007	Solexa-15903	Unrchd Box #3	Solexa-15914	Enrchd Box #7	Solexa-15925	Norm Box #11	3.88
7	E. coli MG1655 WGS	Holly	10/8/2007	Solexa-15904	Unrchd Box #3	Solexa-15915	Enrchd Box #7	Solexa-15926	Norm Box #11	3.70
8	E. coli MG1655 WGS	Holly	10/8/2007	Solexa-15905	Unrchd Box #3	Solexa-15916	Enrchd Box #7	Solexa-15927	Norm Box #11	3.68
9	Myeloma Tumor	Eddie	10/10/2007	Solexa-15906	Unrchd Box #3	Solexa-15917	Enrchd Box #7	Solexa-15928	Norm Box #12	3.72
10	Myeloma Normal	Eddie	10/10/2007	Solexa-15907	Unrchd Box #3	Solexa-15918	Enrchd Box #7	Solexa-15929	Norm Box #12	3.82
11	MRSA WGS	Danielle	10/15/2007	Solexa-15908	Unrchd Box #4	Solexa-15919	Enrchd Box #7	Solexa-15930	Norm Box #12	3.76
12	MRSA WGS	Danielle	10/15/2007	Solexa-15909	Unrchd Box #4	Solexa-15920	Enrchd Box #7	Solexa-15931	Norm Box #12	3.38
13	MRSA WGS	Danielle	10/15/2007	Solexa-15910	Unrchd Box #4	Solexa-15921	Enrchd Box #7	Solexa-15932	Norm Box #12	3.69
14	MRSA WGS	Danielle	10/15/2007	Solexa-15911	Unrchd Box #4	Solexa-15922	Enrchd Box #7	Solexa-15933	Norm Box #12	3.83

- ▶ Implementing a library numbering system allows for easy tracking, even without barcodes
- ▶ Keeping Excel tracking sheets can associate libraries to each other and to project, assisting in downstream troubleshooting
  - ▶ Records who made the libraries, prep date, lot #s, etc.
  - ▶ Identifies where in the freezer the tubes are stored

# Shearing of Genomic DNA

# Shearing of Genomic DNA




# DNA Shearing with Covaris

- ▶ Founded in 1999, Covaris, Inc. developed a line of instruments capable of providing controlled delivery of acoustic energy to closed vessels.
- ▶ Originally, Covaris targeted the chemical industry in order to facilitate mixing, homogenizing, or pulverization of chemical compounds.
- ▶ We were the first to work with Covaris to develop methods to apply their “adaptive focused acoustics” technology to the process of shearing DNA for next generation sequencing applications.
- ▶ Main reason for moving to Covaris was the throughput:
  - ▶ Up to 96 samples can be loaded into the machine at a time in the E210 plate based model.

# DNA Shearing with Covaris

## ► How it works:

- Acoustic ultrasonic waves aimed at glass shearing vessels containing DNA cause “cavitation” within the sample.
- This extremely rapid creation and high-energy collapsing of tiny bubbles leads to the mechanical shearing of DNA.
- We will see a demo in the Lab Section.



E210 Plate Based Model





S2 Single Tube Based Model



# DNA Shearing with Covaris

- ▶ We have performed a multitude of experiments to understand how the parameters can be manipulated to produce DNA in various sizes:

Desired Fragment Size	Volume	Time (T)	Duty Cycle	Intensity	Cycles Per Burst (CPB)	Z Axis (mm)
150bp	100	300 sec	20%	5	200	0mm
200bp	100	90 sec	20%	5	200	0mm
200bp	50	360 sec	5%	5	200	0mm
300-700bp	100	30 sec	20%	5	200	0mm
300-700bp	50	300 sec	1%	5	200	0mm



- ▶ *You will see more examples of these data in the lab section.*


# Choosing a DNA Shearing Method

	DNA in Desired Range	Reproducible	Customizable Settings	No Additional Materials Cost	No Additional Instrument Required	Multiple Samples at a Time	Compatible with automation
<b>Nebulizer</b>	✓	✓	<i>Limited - can vary time and pressure</i>	✓	✓		
<b>Covaris S2</b> Single Tube	✓	✓	✓				
<b>Covaris E210</b> Plate Based	✓	✓	✓			✓	✓

# Reaction Cleanups



# Reaction Cleanups




# Reaction Cleanup Methods

- ▶ Column based cleanups with Qiagen MinElute columns
  - ▶ Reliable and reproducible, ~80% recovery
  - ▶ Easy to perform
  - ▶ Manual - ideal for lower throughput
    - ▶ Not compatible with automation: reasonable set size <12
    - ▶ Tubes more prone to sample swapping than samples in a plate
- ▶ Magnetic bead based cleanup methods
  - ▶ Reliable and reproducible, 80-90% recovery
  - ▶ Easy to perform
  - ▶ Can be done manual or automated
    - ▶ Manual: small sets of 8-12 can be reasonably done by hand in tubes
    - ▶ Automated: up to 96 can be done at a time using a liquid handler
- ▶ **Broad has moved away from column cleanups to magnetic bead based cleanups to accommodate higher throughput**

# SPRI: An Alternative to Column Cleanups


- ▶ **SPRI: Solid Phase Reversible Immobilization** (*AMPure by Agencourt*)
  - SPRI technology is a magnetic bead based DNA clean up method
  - Beads are specific for double stranded DNA
  - Enzymes, excess primer, etc. are washed away & DNA is then eluted
  - Reference: Hawkins TL et al. Nucleic Acids Res. 1994 Oct 25;22(21):4543-4.
- ▶ How the SPRI chemistry works:



Courtesy of Agencourt

# SPRI: An Alternative to Column Cleanups

- ▶ The length of fragments that SPRI binds can be controlled by varying the amount that is added to the reaction
- ▶ The “X” of SPRI refers to the volume to volume ratio of SPRI to DNA
  - ▶ 1X SPRI is a 1:1 vol of SPRI:DNA
- ▶ The **lower** the X of SPRI, the **higher** the size cut off becomes
  - ▶ See titration at right
  - ▶ 0.7X SPRI removes fragments below 150 bp such as adapter dimer
- ▶ As we will see in the size selection section, this property can be used to perform gel free library sizing




# Application of SPRI - High Throughput Libraries

- ▶ To fill our sequencing capacity, we needed to be able to do the following:
  - ▶ Make 96 libraries or more per week
  - ▶ Start with 3 ug of DNA
  - ▶ Average final fragment size ~300 bp but tight sizing NOT required
- ▶ Implementing Covaris shearing and SPRI cleanups made this possible:
  - ▶ Covaris shear of 96 samples
  - ▶ Automated SPRI cleanups on the Velocity 11 Bravo liquid handler
  - ▶ No size selection method
    - ▶ 150 bp Covaris shearing
    - ▶ 0.7X SPRI to remove dimer after adapter ligation




Completed Library Agilent QC




# Size Selection Methods

# Size Selection



# Insert Size vs. Library Fragment Size





# Choosing a Library Fragment Size

- ▶ A library fragment size of ~400 bp is ideal for optimal cluster density as well as for sequencer optics and cluster-finding algorithms
- ▶ However, the sequencer can tolerate fragment sizes up to 800 bp
  - ▶ Fragment size can be tailored to suit the needs of a given project
- ▶ Desired read length determines the minimum insert size
  - ▶ For paired end, the targeted insert size should usually be at least twice the read length to avoid reading into the adapter or sequencing the same stretch
    - For 36 base paired sequencing: at least 72 bases
    - For 76 base paired sequencing: at least 152 bases
- ▶ **Larger fragments → Larger clusters → Lower max density**
- ▶ **Wider range → Wider variety of cluster sizes → Lower max density**


# Choosing a Size Selection Method

- ▶ Option 1: Agarose gels
- ▶ Option 2: Double-sided SPRI
- ▶ Option 3: No size selection

# Option 1: Agarose Gels

- ▶ Gel size selection is necessary to achieve size deviations of  $< 30\%$ .
- ▶ Broad's recommendations for gel conditions:
  - For fragments  $< 400$  bp: 4% NuSieve 3:1 TAE gel, 85 volts, 3 hours, room temp
  - For fragments  $> 400$  bp: 2% SeaPlaque TAE gel, 110 volts, 3 hours,  $4^{\circ}\text{C}$
- ▶ By taking  $\sim 2$  mm gel slices at desired size range, a size distribution of  $\pm 10\%$  or less can be achieved:



We will see a lab demo of a gel cutting tool we have developed to make these tight cuts easier.



Example of a tight cut library targeting fragment size of 490 bp.


# Application of Tight Cuts - 180 bp Libraries

- ▶ We were requested by our Assembly team to make libraries with inserts that were EXACTLY 180 bp with less than a 10% deviation
  - ▶ 2x101 run = virtual 180 bp read
  - ▶ 180 bp virtual read more robust to assemble than 100 bp alone
- ▶ To accomplish this, we take 3 extremely tight gel slices at 246 bp after adapter ligation (180 bp insert + 66 bp adapters)
  - ▶ Use a gel cutting tool that holds 4 razor blades - Will see a demo in lab section
  - ▶ Perform “high diversity” PCR protocol
  - ▶ QC size & distribution by Agilent chip
    - ▶ After PCR, should be 299 bp
- ▶ Select library closest to 299 bp with smallest standard deviation for sequencing



## Option 2: Double-Sided SPRI


- As discussed previously, the size of fragments that SPRI binds can be varied by varying the “X” to SPRI that is added to the reaction containing the DNA:




- We can utilize this property to perform a gel-free and automated size selection that we have termed “Double-sided SPRI”

## Option 2: Double-Sided SPRI

- By implementing a combination of good shearing with SPRI and “reverse” SPRI, one can select a fairly tight size range *with no gel*:




Results:



# Option 3: No Size Selection At All

- ▶ If tight sizing not required for a project, size selection may be unnecessary.
- ▶ Combining a tight 150 bp shear with a 0.7X SPRI clean-up after adapter ligation to remove excess adapter dimer may be suitable.


**POST SHEARING** (150 bp protocol)



End Repair  
A Addition  
Adapter Lig  
NO GEL  
0.7X SPRI

PCR Enrichment  
1 ul DNA  
10 cycles

**POST PCR ENRICHMENT & 0.7X SPRI**




- ▶ This works because of the tightness of our 150 bp shear. If shearing isn't optimized or if a larger size is needed, than double-sided SPRI is required.

# PCR Enrichment



# PCR Enrichment



# Library Complexity

- ▶ Library complexity refers to the number of unique fragments present in a given library
- ▶ Complexity is affected by:
  - ▶ Amount of starting material
  - ▶ Amount of DNA lost during cleanups and size selection
  - ▶ Amount of duplication introduced via PCR
- ▶ For most libraries that only need to be run across a few lanes, the standard protocol provides libraries with ample complexity
- ▶ However, certain projects require very deep coverage from a single sample - i.e. SNP discovery, mammalian assembly, cancer resequencing
- ▶ When dozens of lanes are required, library complexity becomes very important


# Relationship Between PCR and Complexity

## Experiment:

- ▶ Varied input volume, PCR cycle number, & number of reactions pooled
- ▶ Sequenced to determine % duplication & estimated library size
- ▶ Based on data, number of new molecules seen after each successive lane of sequencing was estimated

## Results:

- ▶ Decreasing PCR cycles decreases % duplication
- ▶ Increasing DNA input and pooling multiple PCR reactions decreases % duplication




Starting DNA	Cleanup Method	DNA Input to PCR (ul)	PCR Cycle #	# of PCR Rxn Pooled	% Duplication	Estimated Library Size
3 ug	MinElute	1	10	1	2.79%	247,441,814
3 ug	MinElute	1	14	1	5.51%	195,858,361
3 ug	MinElute	6	10	1	0.83%	1,255,845,288
3 ug	MinElute	1	10	6	0.48%	1,390,926,598

# Application - “High Diversity” Sample Preparation


- ▶ Cancer resequencing projects required 30X unique coverage from tumors and their matched normals
- ▶ We needed a way to reach this goal while minimizing both the amount of genomic DNA used and the amount of sequencing lanes needed
  - ▶ Highest complexity, lowest % duplication
  - ▶ Maintain a very tight size distribution of about 500 bp +/- 10%
- ▶ To accomplish this, we developed a high diversity protocol:
  - ▶ Take **two** extremely tight gel cuts from each sample: 500bp & 520bp
    - ▶ Creates two libraries for each 3 ug input
    - ▶ Double the amount of fragments from the same starting input
  - ▶ Enrich **ENTIRE** sample from both cuts to maximize diversity
    - ▶ 16 x 50 ul PCR reactions with 10 cycles enrichment
    - ▶ Pool 16 reactions for each cut & clean up in same column or SPRI

# Application - “High Diversity” Sample Preparation






# Library Quality Control

# Library Quality Control



# Library Quality Control





- ▶ **Agilent BioAnalyzer 2100**
  - Ensure proper size distribution after both shearing & PCR
  - Failing samples enter rework loop: re-sheared or re-enriched
- ▶ **Pico Quantification Post-Enrichment**
  - Samples <20nM are re-enriched with more input and/or cycles
- ▶ **Gel Images**
  - Look for excessive adapter, DNA in unexpected size range, etc.
  - Pictures taken before and after cutting are saved for downstream troubleshooting






# Library Quality Sequencing Metrics

- ▶ **Always check the sequencing metrics for a run before looking at the LC metrics**

Metric	Definition	How to Interpret...
% Aligned	% of passing filter (PF) reads that aligned to reference sequence.	Low % aligned may indicate sample contamination or swap.
% Adapter 	% of PF reads whose first 16 bases match any part of the Illumina adapter sequence.	>1% adapter indicates inefficient removal of adapter dimer in size selection.
% Chimerism 	% PF of reads that have 2 ends over 100kb apart or on 2 chromosomes.	>1% chimerism indicative of problem in adapter ligation or with genomic DNA prep itself.
% Duplication 	% of PF aligned reads originating from duplicate fragments (i.e. multiple reads with exact same R1 and R2 start sites).	High % duplication indicates a low complexity library, possibly due to low amount of starting material and/or excessive PCR cycles.
Estimated Library Size 	Estimated number of unique molecules in library, calculated using % duplication and reference genome.	Library size depends on starting material. Human WGS libraries should have a size of 1-3 billion.

# Library Quality Sequencing Metrics

Metric	Definition	How to Use Interpret...
Concordance	Concordance of sequencing results with previous genotyping/SNP data.	Can detect sample swap or contamination events.
Insert Size Plots	Show the size distribution of sequenced fragments.	“Shoulders” or spurious peaks indicate problem in size selection.



# Module Summary - Sample Preparation

- ▶ Sample preparation is the most customizable area of the Illumina process - Can adapt the basic protocol to meet the goals of any project
- ▶ Sample tracking is important for any size laboratory to ensure sample integrity
  - ▶ High throughput: barcoded libraries, plates, and LIMS tracking
  - ▶ Low throughput: Excel sheets or databases to record project details & QC information
- ▶ Broad has developed many modifications to the standard protocol to increase throughput and reproducibility, as well as to meet specific project needs:
  - ▶ Optimized shearing with Covaris
  - ▶ Implementation of AMPure SPRI for clean-ups and gel-free size selection
  - ▶ Optimized gel cutting conditions for libraries with <10% size deviation
  - ▶ Development of high diversity protocols for deep coverage libraries
- ▶ Library quality control monitoring and sample tracking are essential to ensure that only high quality libraries make it to the sequencer

# What You Will See in the Lab Sections:

- ▶ In lab sample tracking methods and tools
- ▶ Covaris E210 demo
- ▶ Agarose gel cutting tips & tricks
- ▶ Library QC methods and review of example data
- ▶ SPRI clean-up hands-on demo & automation overview